

## Emergency Response at Critical Control Points

### Point of Use Emergency Kit

- Retain bacteria in water for washing and drinking to aid in infection control
- Produce water that is suitable for cleaning of equipment used in medical procedures and washing of surgeon's hands
- Provide a 90-day FDA 510k-cleared solution to provide time for remediation and install a programmatic filtration solution

Contact us for more information

(201) 343-5202

sales@nephros.com

**An instant, patient-centered solution during compromised water events.**



#### For Emergency Response at Sinks, Showers, or any Point of Use Applications

- FDA-cleared for 90 days
- Contents:
  - (4) S100 Sink Spout Microfilters
  - (4) S100 Shower Handheld Microfilters
  - (4) CPC Sink Adapters (for sinks)

#### Product Number

70-0297F Point of Use Emergency Kit

### Specifications

<b>Max Inlet Pressure</b>	75 psi (5 bar)
<b>Filter Membrane</b>	Versatile® - PES
<b>Material</b>	Polysulfone
<b>Pore Size</b>	0.1 micron
<b>Bacteria Retention</b>	>10 <sup>9</sup> (B. diminuta)
<b>Virus Retention</b>	n/a
<b>Dimensions</b>	Sink: 4" L x 3" D Shower: 4" L x 3" D head + 7" handle
<b>Surface Area</b>	0.3 m <sup>2</sup>
<b>Clean Water Flow Rate</b>	Sink: 1.5 GPM at 50 psi (3.5 bar) Shower: 3.0 GPM at 30 psi (2.0 bar)
<b>Connections</b>	Sink: QC (CPC (PLC)) Shower: Threaded (1/2 NPSM)
<b>Replacement</b>	Up to 3 months (sinks)
<b>Configuration</b>	Single-Stage Filtration

(1) Versatile® is a registered trademark of Medica S.p.A.

NEPHROS.COM

Nephros is committed to improving the human relationship with water through leading, accessible technology. We provide innovative filtration and pathogen detection as part of an integrated approach to water safety that combines science, solutions, and support services. Nephros products serve the needs of customers within the healthcare and commercial markets, offering both proactive and emergency responses for water management.