



Questions or issues?
Let us help!

Call:
(201) 343-5202, Ext. 1

Installation Guide Part #: **70-0285**

70-0286

DSU-H TEMPLATE

Allow 8 inch clearance above and below filter connectors to provide enough space for flex hose connections. Must be printed on 8.5 x 14 paper with no scaling.

MANUFACTURER:

Nephros, Inc.
380 Lackawanna Place,
South Orange, NJ 07079 USA
(201) 343-5202

ASSEMBLED AT:

MEDICA S.p.A.
Via Degli Artigiani, 7
41036 Medolla (MO), Italy

NEPHROS.COM