

Installation Kit – S100 In-Line Shower Head & Arm⁽¹⁾

70-0350

This installation kit is for use with the Nephros S100 In-Line Threaded Microfilter, allowing an infection control filter to be used as a fixture filter in showers with fixed water supply pipes. The arm lifts the filter above the wall outlet to preserve usable vertical shower space. The kit also allows for flexible positioning of the filter, reducing the risk of product breakage.

Kit Contents

- (1) Adjustable Spray Shower Head
- (1) Fixed Shower Adapter Arm

Contact us for more information

(201) 343-5202

sales@nephros.com



⁽¹⁾ Kit components are sourced from third parties and assembled for convenience. Kit components are not manufactured by Nephros.

NEPHROS.COM

Nephros is committed to improving the human relationship with water through leading, accessible technology. We provide innovative water filtration products and services, along with water-quality education, as part of an integrated approach to water safety. Nephros goods serve the needs of customers within medical and commercial markets, offering both proactive and emergency solutions for water management.

60-8056 REV00